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## ACIDITY ON OIL

Turn on the Minifood.

Once warm up is completed, insert the test tube with the reagent in the incubation cell for 3-5 minutes and then press **ACID** to start Low Acidity test. Press again ACID to choose High Acidity test.

Agitate the test tube, inverting it.

Insert the test tube with reagent in the reading cell and press **ENTER**.

Open the test tube.

Important: place the cap with the internal part facing up. Do not touch the internal part with the hands.


Agitate the sample several times to make it homogenous.

Take 2.5  $\mu\text{l}$  of sample for Low Acidity or 1.25  $\mu\text{l}$  for High Acidity: immerse the tip of the pipette into the oil with the piston lowered. **SLOWLY** lift the piston keeping the top of the tip immersed in the sample.

Clean the tip using adsorbent paper as indicated in the photo (keep the pipette horizontal during this operation).

Immerse the tip into the liquid in the test tube and pump the piston up and down 2-3 times by using the thumb.

Close the test tube with its cap and agitate, inverting it.

Insert the test tube in the reading cell and press **ENTER**.

The result is displayed immediately.

## PEROXIDES ON OIL

Turn on the Minifood.

Once warm up is completed, insert the test tube with the reagent in the incubation cell for 3-5 minutes and then press **PEROX** to start Peroxides test.

Agitate the test tube, inverting it.

Insert the test tube with reagent in the reading cell and press **ENTER**.

Open the test tube.

Important: place the cap with the internal part facing up. Do not touch the internal part with the hands.


Agitate the sample several times to make it homogenous.


Take 2.5  $\mu\text{l}$  of the sample: immerse the tip of the pipette into the oil with the piston lowered. SLOWLY lift the piston keeping the top of the tip immersed in the sample.

Clean the tip using adsorbent paper as indicated in the photo (keep the pipette horizontal during this operation).

Immerse the tip into the liquid in the test tube and pump the piston up and down 2-3 times by using the thumb.

Close the test tube with its cap and agitate, inverting it.

Open the test tube and insert one drop of R2 reagent.

Close the test tube again and agitate, inverting it.

Insert the test tube in the reading cell and press **ENTER**.

The result will be displayed after 2 minutes.


---

**CORRECT PROCEDURE****INCORRECT PROCEDURE**

---

The caps has be opened carefully by putting pressure to the right and then to the left. Do not touch the internal parts (soaked) of the cap with the hands.

Do not open the cuvette using force: it will never open.

---

When taking the sample always agitate it, inverting the container with the sample and immerse the tip keeping the piston lowered.

Do not raise the piston fast but take the sample SLOWLY.

---

It is fundamental for making correct analysis that the tip is clean outside: clean it with adsorbent paper as shown.

Do not touch the top of the tip in order not to take the liquid from inside.

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