SLB222 analyser by CDR - Short Guide


Short instructions do not replace the complete user's manual.

Analysis systems based on the SLB222 analyser are:

CDR FoodLab®, CDR WineLab®, CDR BeerLab®, CDR CiderLab, CDR PalmOilTester, CDR OxiTester, CDR KombuchaLab, CDR CoffeeLab®


Rear panel


For the installation of the instrument, proceed as follows:

- A.** Place the **SLB222 analyser** on a flat surface in an area that it is not exposed to mechanical vibrations or extreme temperature fluctuations or excessive light.
- B.** Connect the external power adapter to the port at the rear of the instrument using the appropriate power cord to an outlet between 110V to 230V.
- C.** Press the ON/OFF button to switch ON the analyser. Press and hold the ON/OFF button for six seconds to switch OFF.


D. After powering ON the analyser, the photometer optics will be controlled together with other functional parts and on the screen are shown:

and then

At the end of the heating process, the **Main dialogue** is shown:


The analyser is now ready to be used.

How to check if the analyser has the latest software version.


We suggest to carry out the check after the first switch-on.

1. Connect the instrument to your PC via the USB-type B cable

2. Find the CDRDRIVE in your PC and open it

3. Double click on the CONNECT link you'll find in the CDRDRIVE and follow the instructions

4. The resources.cdrfoodlab.com website will open and your analyser will be recognized. If it is suggested on the webpage, update the analyser following the instructions.

CDR FoodLab TECHNICAL SUPPORT

Dashboard | Analysis/MSDS | Manuals | Software Updates | CDR Data Manager | Warranty | Applications

Warranty Information
Please, upload invoice of your instrument to extend the warranty period to 36 months.

Firmware Information
Your instrument is out of date. Update it!

CDR FoodLab

Series	222M7
SN	1238
Version	10.19
Register by	Pierpaolo Sartori
Company	Lasteria Montello

Last connection 11/09/2024 11:58:01

MSDS - Quick access to the security information
Search the security information using the code on the package of your reagents kit.

Analysis Certificate
Request the certificate of analysis for the reagents kit:

REF	REF
LOT	LOT
Email	

Analysis / MSDS
Discover the technical specifications and resources related to the analyses of your instrument.

Manuals
Download quick start guides and manuals to discover all the instrument's features. Here you can find all the documents you need.

CDR SLB222 Installing printing paper

Pull up the green lever on the printer to release the black roller and insert the paper under the printer roller.

Pull down the green lever to secure the printing roll and turn the green wheel to adjust the paper.

FOOD TEST PERFORMANCE EQUIPMENT

Rev. 1.1 of 18.03.2015

WARNING!

Use the device by strictly following the rules and requirements listed below

- Before using the device, all operators should be familiar with the instructions for use that can be downloaded from CDR website.
- Before the first operation, it is recommended that the packaging and the device be visually checked to verify that it has not been damaged during transport. Return the device to the Maintenance department for a functional inspection if damage is reported.
- The device is designed to perform food tests and can be used with cuvettes, code CDR 230012, only.
- The device cannot be used for tests on biological matrices that differ from those described in methods and for which the device has been specifically calibrated during manufacturing.
- The device cannot be used with assay samples that have not been specifically released by CDR.
- The device cannot be used to perform continuous diagnostic tests.
- Do not place reflecting objects inside the wells, nor get close to the well opening with your eyes to avoid exposure to light radiations, which are considered, on a prudential basis, as potentially hazardous to the sight.
- The manufacturer recommends to preserve the original box in the event it is necessary to transport the device and install it in another location. In terms of disposal, cardboard is not regarded hazardous for people, animals or materials. Before disposing of the package, see Directive 94/62/EC on packaging and packaging waste.

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